

MEDICAL & HOSPITAL



Under directed temperature



Temperatures control for blood pockets



The blood is sensitive to the variations of temperature. In case of break of the cold chain, it becomes dangerous for health. "GIOVATEMP" is a proof of the preservation temperature of your products. If the maximum temperature allowed is exceeded, a red cursor progresses in the middle of the indicator and will indicate you the approximate overrun time of your products. This will allow you to know simply if your blood pockets can be transfused or not.

Temperatures control for transplants



Auto-adhesive capsule

For a graft or an organ transport from one institution to another, « AURETEMP » allows you to be signalled of any incident of temperature. If the maximum temperature allowed is exceeded, the indicator will trend to red on its whole surface for a preset time and will indicate you if your organs can be implanted or not.

C.C.F

TECHNOLOGIES

www.ccf-technologies.com